



# SINGLE-USE INSTRUMENTS CATALOGUE

[www.mddiglobal.com](http://www.mddiglobal.com)

**mddi** |   
STERILE | SINGLE USE

# CONTENTS

OUR INSTRUMENTS	3
.....	
INDIVIDUAL DIAGNOSTIC INSTRUMENTS	4
.....	
INDIVIDUAL HYGIENIST INSTRUMENTS	7
.....	
INDIVIDUAL RESTORATIVE INSTRUMENTS	11
.....	
DIAGNOSTIC KITS - CPITN-C PROBE	15
.....	
DIAGNOSTIC KITS - PERIO PROBE/SICKLE	18
.....	
RESTORATIVE KITS	20
.....	





## Our Instruments are:

- Quick and easy to set up and use
- Prewrapped in a range of convenient kits and instruments
- Sterile at the point of use
- Convenient to dispose of in sharps container
- Reducing the risk of sharps injuries
- Sharp where needed, ensuring the highest quality of treatment for every procedure
- Implemented to significantly reduce the risk of cross-contamination
- GS1 compliant with removable tracking labels
- CE certified and TGA registered
- 5 year expiry from point of manufacture

## MDDI Can Increase Revenue by:

- Minimise decontamination cycles per shift - reduces expenditure on consumables, test kits and utilities
- Putting the RDN back chairside to assist the dentist and patient
- Increasing the number of patients seen in each clinical session
- Improving clinical efficiencies allowing for longer treatment plans and procedures
- Improving the patient journey





Our *diagnostic instruments* are high quality, precision made instruments. The probes and London College Tweezers are manufactured from stainless steel and all of our handles are recyclable ABS plastic.

View our [individual diagnostic range](#):



# STERILE, SINGLE-USE INDIVIDUAL INSTRUMENTS

## Front Face Mirror with Percussion Tip

Sterile Mouth Mirror - Front Surface Plane 4

SI053ST-1 - 180 per box



## London College Tweezers


SM182ST-1 - 20 per box

SI182ST-1 - 200 per box




## Tweezers

SI081ST-1 - 200 per box



## No. 8 Standard Probe/CPITN-C Probe 3.5/5.5/8.5/11.5

SM084ST-1 - 20 per box



## Perio Probe/Sickle

SM089ST-1 - 20 per box



BEST SELLER

## Williams Probe 10mm/CPITN-C Probe 3.5/5.5/8.5/11.5

SM090ST-1 - 20 per box





Our *hygienist instruments* feature a substantial, ergonomic handle (11mm diameter) with ridges and grooves for improved comfort and accentuated grip.

View our [individual hygienist range](#):

Moulded around a strong, stainless-steel stem which provides durability and weight to ensure the instrument is fit for purpose. Sharp, precision manufactured tips ensure consistent quality for every patient.



## Sickle Scaler 204S

SM202ST-1 - 20 per box

Targets posterior regions of mouth

Removes supragingival interproximal calculus



## Hygienist Scaler H6/H7

SM203ST-1 - 20 per box

Targets inter-proximal anterior and premolar regions



## Gracey Curette 11/12

SM211ST-1 - 20 per box

Targets mesial surfaces of posterior teeth

Removes subgingival calculus



## Gracey Curette 17/18

SM217ST-1 - 20 per box

Targets periodontal pockets and distal surfaces of posterior teeth.



## Gracey Curette 13/14

SM213ST-1 - 20 per box

Target's distal surfaces and pockets of upper/lower premolars and molars

Effective in subgingival debridement





## Universal Curette 2R/2L

SM222ST-1 - 20 per box

Removes moderate to heavy calculus from supra- and subgingival surfaces

## Universal Curette 4R/4L


SM224ST-1 - 20 per box

Targets anterior and posterior teeth

Effective in supra- and subgingival scaling









Our *restorative instrument* handles are moulded around the stainless steel core, ensuring the instruments display strength and rigidity when required.

View our [individual restorative range](#):



### Flat Plastic DE

SM154ST-1 - 20 per box




### Medium Spoon Excavator DE (1.8mm)

SM158ST-1 - 20 per box



### Ball Burnisher DE

SM156ST-1 - 20 per box





## Amalgam Plugger Large/Small

SM161ST-1 - 20 per box



## Wards Carver No.2

SM162ST-1 - 20 per box



## Spatula DE

SM165ST-1 - 20 per box



## Ladmore Plugger DE

SM157ST-1 - 20 per box



## Machtou Plugger 0.5mm/0.6mm

SM163ST-1 - 20 per box



## Machtou Plugger 0.8mm/1.0mm

SM164ST-1 - 20 per box



# STERILE, SINGLE-USE DIAGNOSTIC KITS



Our MDDI diagnostic kits were adopted into the following sectors more than 15 years ago. We have continued to expand our range due to the ongoing popularity of our instruments.


School  
Clinics

Government  
Hospitals

Flying Doctors  
Service

Mobile  
Dental Units





## Periodontal 2 Piece Examination Kit with CPITN-C Probe

SM952ST-1 - 25 kits per box  
SI952ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip  
1 x No.8 Standard Probe/CPITN-C Probe

BEST SELLER

## Periodontal 3 Piece Examination Kit with CPITN-C Probe

SM951ST-1 - 25 kits per box  
SI951ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip  
1 x No.8 Standard Probe/CPITN-C Probe  
1 x Tweezers



## Premium Periodontal Examination Kit with CPITN-C Probe

SM632ST-1 - 25 kits per box  
SI632ST-1 - 100 kits per box

1 x Front Face Mirror with Percussion Tip  
1 x No.8 Standard Probe/CPITN-C Probe  
1 x London College Tweezers





## Orthodontic Three Piece Examination Kit with CPITN-C Probe

SI634ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
- 1 x Orthodontic Ruler (150mm)



## Suture Removal Kit

SI770ST-1 - 25 kits per box



- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
- 1 x Iris Scissors Curved 11.5cm
- 1 x Iris Toothed Forceps 10cm

- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 2 x Non-Woven Swab 4ply (5cm x 5cm)
- 1 x Tray (19.5cm x 12.5cm x 2cm)



## Complete Periodontal Examination Kit with CPITN-C Probe

SI950ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x No.8 Standard Probe/CPITN-C Probe
- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)



## PERIO PROBE/SICKLE EXAMINATION KITS



### Periodontal 2 Piece Examination Kit

SM552ST-1 - 25 kits per box  
SI552ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle

### Periodontal 3 Piece Examination Kit

SM551ST-1 - 25 kits per box  
SI551ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers





## Premium Periodontal 3 Piece Examination Kit

SM631ST-1 - 25 kits per box  
 SI631ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x London College Tweezers

## Orthodontic Three Piece Examination Kit

SI633ST-1 - 100 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Orthodontic Ruler (150mm)



## Complete Periodontal Examination Kit

SI550ST-1 - 25 kits per box

- 1 x Front Face Mirror with Percussion Tip
- 1 x Perio Probe/Sickle
- 1 x Tweezers
- 2 x Type IIR Face Mask with Ear Loops
- 1 x Patient Bib with Ties
- 2 x Cotton Roll No.2
- 1 x Tissue Wipe
- 1 x Tray (19.5cm x 12.5cm x 2cm)





# STERILE, SINGLE-USE RESTORATIVE KITS



Labour costs make up a significant portion of a clinic's annual budget with sterilisation of equipment a major contributor to the daily workload.

**This is where our kits can make a difference.**

Save time in the steri room

+ Save money on equipment and supplies

+ Improve infection control

Our comprehensive range of restorative kits are convenient, sharp and new every time.



## 2 Piece Restorative Mixing Kit

SI622ST-1  
150 kits per box

1 x Spatula DE  
1 x Dycal Applicator



## 2 Piece Composite Restorative Kit

SI620ST-1  
150 kits per box

1 x Flat Plastic DE  
1 x Ball Burnisher DE



## 3 Piece Composite Restorative Kit

SI630ST-1 - 100 kits per box

1 x Flat Plastic DE  
1 x Ball Burnisher DE  
1 x Medium Spoon Excavator DE (1.8mm)



BEST SELLER



## 4 Piece Composite Restorative Kit with London College Tweezers

SI640ST-1 - 80 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x London College Tweezers



## 5 Piece Amalgam Restorative Kit

SI650ST-1 - 60 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x Carver DE
- 1 x Amalgam Plugger Large/Small



## 6 Piece Composite Restorative Kit

SI660ST-1 - 50 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Small Spoon Excavator DE (1.2mm)
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x PFI49 Burnisher
- 1 x 1/2 Hollenback DE



## 6 Piece Comprehensive Restorative Kit

SI661ST-1 - 50 kits per box

- 1 x Flat Plastic DE
- 1 x Ball Burnisher DE
- 1 x Medium Spoon Excavator DE (1.8mm)
- 1 x Carver DE
- 1 x Ladmore Plugger DE
- 1 x Amalgam Plugger Large/Small



MDDI instruments are available to purchase through Europe's leading distributors.

Visit our [distributor page](#) to find out more.



 HENRY SCHEIN®



We are constantly looking for new ways to collaborate, please do not hesitate to get in touch.

**CONTACT MDDI UK:**

P: +44 (0)1535 657 209 | E: [salesuk@mddiglobal.com](mailto:salesuk@mddiglobal.com)

1 Marrtree Business Park, Ryefield Way, Silsden, West Yorkshire BD20 0EF UK

**MDDI - HEAD OFFICE**

E: [sales@mddiglobal.com](mailto:sales@mddiglobal.com) | P: +61 (0)7 5453 4433

Unit 1, 24 Hoopers Road, Kunda Park Queensland 4556 AUSTRALIA



[www.mddiglobal.com](http://www.mddiglobal.com)